

# **Omni Bladder Management System**

## **MFG Part: CCU-10**

### **Product Name: Waist Collection Chamber Unit**

#### **ProRen – Prosthetic Renal Management System**

### **Drainable Collection Chamber Unit**

- Capacity of 1000mL
- Made from durable, medical-grade urethane
- Latex-free
- Patented Control Unit mounting clip
- Air filtering, hydro block vents
- Twist valve
- Anti-reflux valve
- Includes quick disconnect belt



The Waist Collection Chamber Unit is used with the System Control Unit and ProRen CUP (with Harness /Brace) or the ProRen FLO to provide non-invasive, ambulatory and non-ambulatory bladder relief that keeps the patient dry.

**When ordering Collection Chamber, specify “Waist Chamber”.**

- Automatically collects urine as the Control Unit pumps it away from the ProRen CUP or ProRen FLO.
- Weight : 1.9 oz
- Dimensions: 12” x 7.5”

These products are protected by one or more of the following United States patents:  
6,706,027 / 6,918,899 / 7,131,964 / 7,135,012 / 7,141,043 / 7,335,189 / 7,866,942 or other patents pending

**Omni Medical Systems**, an Omni Measurement Systems, Inc. company  
P: 802-497-2253 F: 802-497-3601 Email: [Sales@OmniMedicalSys.com](mailto:Sales@OmniMedicalSys.com) [www.OmniBladder.com](http://www.OmniBladder.com)

