

# ***Omni Bladder Management System***

***MFG Part: CCU-10***

***Product Name: Collection Chamber Unit***

***ProRen – Prosthetic Renal Management System***

## **Drainable Collection Chamber Unit**

- Capacity of 600mL, 1000mL, or 2000mL
- Made from durable, medical-grade urethane
- Latex-free
- Patented Control Unit mounting clip
- Air filtering hydro block vents
- T-valve for easy draining
- Anti-reflux valve



The Collection Chamber can be worn sitting, standing, walking or semi-reclined.

The Collection Chamber Unit is used with the System Control Unit and ProRen CUP (w/Harness/Brace) or ProRen FLO to provide non-invasive, ambulatory and non-ambulatory bladder relief that keeps the patient dry.

**When ordering, please specify capacity.**

- Automatically collects urine as the Control Unit pumps it away from the Male Cup.
- Weight : 1.9 oz
- Dimensions: 10" x 13"

These products are protected by one or more of the following United States patents:

6,706,027 / 6,918,899 / 7,131,964 / 7,135,012 / 7,141,043 / 7,335,189 / 7,866,942 or other patents pending

**Omni Medical Systems**, an Omni Measurement Systems, Inc. company

P: 802-497-2253 F: 802-497-3601 Email: [Sales@OmniMedicalSys.com](mailto:Sales@OmniMedicalSys.com) [www.OmniBladder.com](http://www.OmniBladder.com)



600mL Chamber

