

***ProRen Bladder Management System***  
***MFG Part: SCU-RBP-B***  
***Product Name: System Collection Unit (ProRen FLO)***

***ProRen – Prosthetic Renal Management System***

**System Collection Unit (ProRen CUP)**

- Urine sensor for automatic operation
- Worn discreetly
- Non invasive urethane cup
- Soft foam ring comfortably fits to body
- Fits like a flexible male athletic cup
- Cup stays comfortably in place and ready to use for ambulatory or wheelchair users



**(When ordering this item Specify 12", 18", 24" or Custom hose length.)**

The ProRen CUP is used with the System Control Unit , Harness/Brace, and Collection Chamber to provide comfortable, non-invasive bladder relief that keeps the patient dry.

- Automatically senses the presence of moisture and activates the system
- Flow rate 1.3LPM +/- .2LPM
- Attaches to Collection Bag for fully-automatic use
- Weight : 4.8 oz
- Approximate size similar to Athletic Cup
- Reusable with cleaning

These products are protected by one or more of the following United States patents:

6,706,027 / 6,918,899 / 7,131,964 / 7,135,012 / 7,141,043 / 7,335,189 / 7,866,942 or other patents pending

**Omni Medical Systems**, an Omni Measurement Systems, Inc. company  
808 Hercules Drive, Colchester, VT 05446 802-497-2253 sales@omnimedicalsys.com

