

# ***ProRen Bladder Management System***

***MFG Part: SCU-RBP-VA***

***Product Name: System Collection Unit (ProRen FLO)***

## **System Collection Unit (ProRen FLO)**

- Hydroblock air vents allow condom to fully drain and keep skin dry
- Non invasive silicone sheath
- Hydrocolloid adhesive
- Significantly reduces adverse events
- Prevents accidental detachment
- Latex and Vinyl free
- FLO Collection Unit stays comfortably in place for ambulatory and non ambulatory users
- Quick Disconnect design allows for intermittent catheter use without removing condom.



**(When ordering this item specify STYLE and DIAMETER.)**

The ProRen FLO is used with the System Control Unit , System Chamber Unit, and Collection Chamber Unit to provide reliable, comfortable and non-invasive bladder relief that keeps the patient dry 24 hours a day.

- Each FLO Collection Chamber can be worn for up to 48 continuous hours
- Automatically senses the presence of moisture and activates the system
- Automatic positive pump flow rate 1.3LPM +/- .2LPM prevent condom ballooning
- Attaches to Collection Bag and Control Unit for fully-automatic use
- The hose and sensor (System Chamber Unit) quick disconnect is reusable with cleaning

These products are protected by one or more of the following United States patents:

6,706,027 / 6,918,899 / 7,131,964 / 7,135,012 / 7,141,043 / 7,335,189 / 7,866,942 or other patents pending

**Omni Medical Systems**, an Omni Measurement Systems, Inc. company

808 Hercules Drive, Colchester, VT 05446 802-497-2253 sales@omnimedicalsys.com



# ProRen Bladder Management System

MFG Part: **SCU-RBP-VA**

Product Name: **System Collection Unit (ProRen FLO)**

## System Collection Unit (ProRen FLO)

- FLO System hydroblock air vents allow condom to fully drain and keep skin dry
- FLO System holds the end of the condom catheter open, preventing it from collapsing



The ProRen FLO uses a silicone condom catheter with a moisture-wicking hydrocolloid adhesive, combined with the FLO System to ensure superior performance.

The ProRen FLO's breathable silicone design works well for a wide range of patients.



### Sizing

Small	25mm
Medium	29mm
Intermediate	32mm
Large	36mm
X-Large	41mm

- Non invasive silicone sheath
- Hydrocolloid adhesive



The **Standard** offers a built-in adhesive band in a length and placement found with traditional male external catheters.



The **Short** provides forward placement adhesive technology with a 1.5 inch sheath, perfect for men who have a short or retracted penis.



The **Full** provides 70% more adhesive area than traditional designs. This design helps prevent accidental catheter detachment.



The **Wrap** is an all-silicone, non-adhesive catheter ideal for patients that cannot use catheters with adhesives. Each box comes with three adjustable and reusable foam straps.

**(When ordering this item specify STYLE and DIAMETER.)**

**Omni Medical Systems**, an Omni Measurement Systems, Inc. company

808 Hercules Drive, Colchester, VT 05446 802-497-2253 sales@omnimedicalsys.com

